


SAFE, TIME-SAVING DISPOSABLES FOR ECT



THYMAPAD™ Adherent Stimulus Electrodes (US pat. 4,870,969)

Thymapads™ are fast and easy to apply and remain exactly where you place them. They are so much faster and easier to use than old-fashioned metal disks, rubber headstraps, and jelly, that we have sold over four million of them so far.

Thymapads™ have no exposed metal surfaces to cause accidental shocks, and there's no mess to clean up afterwards, no headstraps or electrodes to wash and dry, nothing to sterilize, no sticky hands to wash—just remove and discard.

#EPAD
MORE THAN
4.0 MILLION SOLD

The patented *Thymapad™* gel surface provides lower impedances and flexibly conforms to the curved surface of the head wherever positioned.

VENTIL-A™ Mouth Protector (US pat. 6,039,046)

The *Ventil-A™*'s thick 100% closed-cell foam construction protects all the teeth, unlike the skimpy cardboard-and-foam product that costs almost 50% more and covers only a few molars.

The *Ventil-A™* easily fits under any anesthesia mask; its patented design features a non-collapsible air channel for free flow of oxygen. Its dimensions are based on measurements of dozens of dental impressions of men and women of all ages, so it fits >98% of adults.



#VENT
MORE THAN
575,000 SOLD

The one-piece *Ventil-A™* design provides dimensional stability; its looped end facilitates fast, easy insertion and removal.

Both *Thymapad™* and *Ventil-A™* single-use ECT aids save the time and expense of washing and sterilization, and eliminate the risk of cross-infection that occurs with re-usable products.

Special Introductory Combination Price for one case of each when ordered with a new Thymatron or under an existing warranty.

Available Exclusively from Somatics: The Place For All Your ECT needs

Call 1(800) 642-6761 Today For a Free Trial Supply



SOMATICS, LLC
910 Sherwood Drive, Unit 23
Lake Bluff, IL 60044

USA & Canada toll-free 1-800-642-6761
(847) 234-6761 www.thymatron.com
FAX (847) 234-6763 sales@thymatron.com

This material was copied at the NLM and may be Subject US Copyright Laws

NATIONAL LIBRARY OF MEDICINE



NLM 03070722 6

Automated.

Easy.

Fast.




Thymatron® System IV

How Smart ECT Just Became Easier to Give

Some feel ECT instruments are too complex. That's why we ship each Thymatron® System IV with a ready-to-use highly effective treatment program enabled so you can just:

Set Dose, Test Impedance, and Treat

This optimal 0.5 msec stimulus program *automatically* chooses the best frequency and duration at every dose. To make it even easier, now you can just press and turn the % Energy dial to switch instantly to the *Old Faithful* stimulus program of the Thymatron® DGx, a 0.25 msec ultrabrief program, or any other stimulus program you prefer, including one you have created and saved.



All the other advanced features can be used or not as you prefer, including unique automatic treatment quality measurements, and EEG, ECG & EMG printouts.

Best of all, the Thymatron® System IV provides the most effective and efficient ECT stimuli available.




SOMATICS, LLC
910 Sherwood Drive, Unit 23
Lake Bluff, IL 60044

USA & Canada toll-free 1-800-642-6761
(847) 234-6761
FAX (847) 234-6763
sales@thymatron.com

www.thymatron.com
This material was copied
at the NLM and may be
Subject US Copyright Laws

Automated.

NATIONAL LIBRARY OF MEDICINE



NLM 04236988 0

**NEW
REMOTE**



**TREAT
HANDLE!**

Easy.

Fast.




Thymatron® System IV

Smart ECT Now Easier than Ever to Give

Do you feel ECT machines are too complex? We do, which is why we ship each Thymatron® System IV preinstalled with a ready-to-use exclusive Optimal Stimulus Program that prevents ineffective treatments by automatically selecting the most effective combination of frequency and duration. The preferred 0.5 ms pulsewidth is uniquely available at every dose, including the maximum dose (or, you can instantly select the 0.3 ms ultra-brief stimulus delivered via the same program). Now you can just:

Set Dose, Test Impedance, and Treat

We've also made it easier to rapidly switch to other stimulus programs—including ones you have created and saved, or any other you prefer—by just pressing and turning the single % Energy dial, preventing the accidental mis-settings that can occur with 4 different dials. And, we have responded to your requests for a Remote Treat Handle by introducing the shielded model shown above that shortens setup time and allows you to deliver the stimulus with a single thumb press.



Best of all, the 900 mA stimulus of the Thymatron® System IV is 50% more effective than 800 mA at the same charge.

DFR302 PR 305
01-20-11 13270 XI ECT Group



This material was copied at the NLM and may be



SOMATICS, LLC
910 Sherwood Drive, Unit 23
Lake Bluff, IL 60044

USA & Canada toll-free 1-800-642-67
(847) 234-67
FAX (847) 234-67
sales@thymatron.c

Automated.

NATIONAL LIBRARY OF MEDICINE



NLM 04405800 3

**NEW
REMOTE**



**TREAT
HANDLE!**

Easy.

Fast.



Thymatron® System IV

Smart ECT Now Easier Than Ever to Give

Do you feel ECT machines are too complex? We do, which is why we ship each Thymatron® System IV preinstalled with a ready-to-use exclusive Optimal Stimulus Program that prevents ineffective treatments by automatically selecting the most effective combination of frequency and duration. The preferred 0.5 ms pulsewidth is uniquely available at every dose, including the maximum dose (or, you can instantly select the 0.3 ms ultra-brief stimulus delivered via the same program). Now you can just:

Set Dose, Test Impedance, and Treat

We've also made it easier to rapidly switch to other stimulus programs—including ones you have created and saved, or any other you prefer—by just pressing and turning the single % Energy dial, preventing the accidental mis-settings that can occur with 4 different dials. And, we have responded to your requests for a Remote Treat Handle by introducing the shielded model shown above that shortens setup time and allows you to deliver the stimulus with a single thumb press.



*Best of all, the 900 mA stimulus of the Thyn
mA at the same charge.*

DFR352 PR 245
04-14-11 13270 XI Group

 **SOMATICS, LLC**
910 Sherwood Drive, Unit 23
Lake Bluff, IL 60044



PROPERTY OF THE
**NATIONAL
LIBRARY OF
MEDICINE**

Free Trial!


www.thymatron.com

This material was copied at the NLM and may be


42-674
(847) 234-674
FAX (847) 234-674
sales@thymatron.co

Somatics: Provider of All Your Advanced ECT Needs for Over a Quarter Century


Ventil-A™ Disposable Oral Protectors



Thymapad™ Disposable Stimulus Electrodes



MouthGuards™




Recording Electrodes



MicroStim Muscle Relaxation Tester




EctoBrain™ ECT Stimulus Analyzer




Remote Treat Handle



The Latest on ECT



Fan-Fold Recording Paper



*Available Exclusively from Somatics:
The Place For All Your ECT needs*



SOMATICS, LLC
910 Sherwood Drive, Unit 23
Lake Bluff, IL 60044

USA & Canada toll-free 1-800-642-6761

(847) 234-6761

www.thymatron.com

FAX (847) 234-6763

sales@thymatron.com

This material was copied at the NLM and may be Subject US Copyright Laws

Automated.

Easy.

Fast.

Thymatron® System IV

Set Dose, Test Impedance, and Treat

NATIONAL LIBRARY OF MEDICINE



NLM 05274265 8



Do you feel ECT machines are too complex? We do, which is why we ship each Thymatron® System IV preinstalled with a ready-to-use exclusive Optimal Stimulus Program that prevents ineffective treatments by automatically selecting the most effective combination of frequency and duration. The preferred 0.5 ms pulsewidth is uniquely available at every dose, including the maximum dose (or, you can instantly select the 0.3 ms ultra-brief stimulus delivered via the same program). Remote treat handle now available.

RAPID TEST OF ECT DEVICE INTEGRITY WITH ECTOBRAIN™

Is your Thymatron® delivering the dose you set? Check it in seconds with the patented EctoBrain™ Analyzer, only from Somatics. With a single button press it performs the same stimulus check used by professional engineers, avoiding the delays and costs that accompany maintenance checks. For complete confidence you can have the ECT nurse perform this check every morning before ECT. Also featured is a separate *Patient Simulator Mode* for training purposes that generates EEG, EMG, and ECG signals from actual patient recordings for demonstrating and testing up to 4 monitor-printer channels of high or low intensity response.



Special Price When Ordered Together With a Thymatron® System IV



SOMATICS, LLC
910 Sherwood Drive, Unit 23
Lake Bluff, IL 60044

**Call For A
Free Trial!**

This material was copied
at the NLM and may be
www.thymatron.com

USA & Canada toll-free 1-800-642-676
(847) 234-676
FAX (847) 234-676
sales@thymatron.com